



Choice of Bearing Surface for Total Hip Replacement Affects Need for Repeat Surgery: A Canadian Perspective

Nicole de Guia, Manager, Joint Replacement & MS Registries
CAHSPR Conference 2014, Toronto, Ontario



Canadian Institute
for Health Information

Institut canadien
d'information sur la santé

Authors and Disclosures

- No disclosures
- CIHI authors: Shirley Chen, Nicole de Guia, Shirley Chen, Vivian Poon, Michael Turner Greg Webster
- Clinical Co-Author: Dr. Eric Bohm, Concordia Hospital, Winnipeg & Chair, CJRR Advisory Committee
- Acknowledgement: Stephanie Cheng, CJRR Sr. Analyst
- Contact: cjrr@cihi.ca

Study Background


- In 2012-2013, there were 47,137 acute care hospitalizations in Canada for all hip replacements (5-year increase: 16.5%). Of these, 9.4% were revision procedures, which has remained stable over the years.
- Joint revision surgeries are more complex than primary procedures and place a heavier burden on the health care system.
- Hospital-related expenditure estimates: \$12,802 for typical hip revision compared to \$8,855 for a primary unilateral THA.*

* Based on CIHI's Patient Cost Estimator. Excludes other expenditures such as physician payments and rehabilitation.

Clinical Introduction

- A THA involves replacing the diseased or damaged hip joint with an artificial joint.
- The femoral head and the acetabular liner (and/or cup) make up the two articulating surfaces of the joint implant (bearing surfaces).
- Bearing surfaces can be made of a variety of materials: metal, ceramic and polyethylene (plastic).
- Surgeons carefully choose the types of bearing surface, considering: patient age, sex and level of physical activity, as well as previous clinical experience.

Example: metal ball/stem



Bearing Surface Types

- Many published studies have found associations between THA bearing surface and increased risk of revision:
 - Breakdown of polyethylene acetabular cup leading to the loosening of the implant.
 - Metallic orthopedic implants generating metal debris from wear and corrosion over time, detectable in the blood, tissue and urine.
- Other international registries (Australia, England and Wales) have conducted similar analyses.
- Both the US FDA and Health Canada have also provided alerts regarding the safety of metal-on-metal hip replacements, leading to increased public awareness.

Study Methodology


- Data sources:
 - Canadian Joint Replacement Registry (CJRR), CIHI
 - Pan-Canadian database that collects data on hip and knee replacement procedures performed across the country. CJRR coverage during study period ranged from 37% to 41%, excluding Quebec.
 - Used to identify primary study cohort and bearing surface material used.
 - Discharge Abstract Database (DAD), CIHI
 - Captures administrative, clinical and demographic information on hospital discharges. Select jurisdictions also use the DAD to capture day surgery.
 - Used to identify revision surgeries among study cohort.
- THA primary cohort linked to DAD to identify any revision surgeries within five years of the primary surgery.
- Covariates: patient sex, age group, bearing surface, procedure year, geographic region and presence of comorbidities at time of primary procedure.

Results – THA Study Cohort

- Cohort: 56,942 primary THAs performed across Canada* between 2003–2004 and 2010–2011.
- Cohort attributes:
 - Most common type of bearing surface material used was metal-on-polyethylene (73%).
 - Metal-on-metal hip replacements (total resurfacing, large-diameter modular or stemmed) made up 9% of all THAs.
 - Patients receiving metal-on-metal THAs were more likely to be younger and male.

* Study did not include procedures from Northwest Territories, Nunavut, Prince Edward Island and Quebec.

Figure 1: Distribution of Bearing Surface Group, THA, 2003–2004 to 2010–2011



N=56,942

Notes:


* Other includes bearing surfaces that were reported as “other”, as well as materials each reported <5% (e.g. ceramic-on-metal and metal-on-ceramic).

MoP: metal-on-polyethylene; MoM: metal-on-metal

Sources:

Canadian Joint Replacement Registry and Discharge Abstract Database, 2003–2004 to 2010–2011, CIHI

Figure 2: Distribution of Bearing Surface Group Over Time, THA, 2003–2004 to 2010–2011




Notes:

* Metal-on-metal include: total resurfacing, large-diameter modular and stemmed subgroups. Metal-on-polyethylene include: standard and cross-linked polyethylene subgroups.

† Other includes bearing surfaces that were reported as “other”, as well as materials each reported <5% (e.g. ceramic-on-metal and metal-on-ceramic).

Sources: Canadian Joint Replacement Registry and Discharge Abstract Database, 2003–2004 to 2010–2011, CIHI

Figure 3: Distribution of Bearing Surface Group by THA Patient Age Group, 2003–2004 to 2010–2011



Notes:

* Metal-on-metal include: total resurfacing, large-diameter modular and stemmed subgroups. Metal-on-polyethylene include: standard and cross-linked polyethylene subgroups.


† Other includes bearing surfaces that were reported as “other”, as well as smaller groups (for example, ceramic-on-metal and metal-on-ceramic).

Excludes 12 cases of unknown patient age.

Results – Early Revisions (Within 5 Years)

- A total of 1,438 first-revision procedures were identified.
- In the years following primary THA, most bearing surfaces had similar cumulative revision rates up until the two-year mark, after which a distinct pattern emerged.
- By the five-year mark, large-diameter metal-on-metal THA had the highest cumulative percentage revision among all groups (5.9%).
- Revisions within one year: primarily issues not related to the mechanics of the implant but rather to medical complications (such as infection, traumatic injury and arthrosis).

Figure 4: Cumulative Percentage Revision by Bearing Surface Group



Notes:

* Other includes bearing surfaces that were reported as “other”, as well as smaller groups (for example, ceramic-on-metal and metal-on-ceramic).

MoP: metal-on-polyethylene

MoM: metal-on-metal

Results – Factors Affecting Risk of Revision

- Factors showing a significantly increased risk ($p < 0.05$), after adjustment:
 - Patients under age 55 had a 1.22 times increased risk compared to those age 65–74.
 - Large-diameter modular metal-on-metal THA had a 1.61 times increased risk compared to metal-on-cross-linked-polyethylene THA.
 - Procedures from Eastern Canada had a 1.35 times increased risk compared to those from Central Canada.
 - Patients with at least one comorbid condition had a 1.47 times increased risk compared to those with none.

Table 1: Cox Proportional Hazards Model for Revision within Five Years After Total Hip Arthroplasty



Risk Factor		Hazard Ratio	95% Confidence Interval	p-value
Age Group [reference: 65–74]				
	< 55	1.22	1.02 – 1.45	0.03*
	55–64	1.09	0.94 – 1.27	0.24
	75 +	1.03	0.90 – 1.18	0.67
Bearing Surface Group [reference: MoP Cross-Linked]				
	MoP: Standard	1.02	0.88 – 1.19	0.79
	MoM: Total Resurfacing	1.02	0.76 – 1.35	0.89
	MoM: Large-Diameter Modular	1.61	1.23 – 2.07	< 0.01*
	MoM: Stemmed	0.89	0.61 – 1.27	0.55
	Ceramic-on-Ceramic	1.05	0.86 – 1.29	0.61
	Ceramic-on-Polyethylene	0.96	0.72 – 1.25	0.74
	Other‡	0.95	0.74 – 1.22	0.74

Table 1: Cox Proportional Hazards Model for Revision within Five Years After Total Hip Arthroplasty (cont'd)



Risk Factor	Hazard Ratio	95% Confidence Interval	p-value
Geographic Region [§] [Reference: Central (Man. Ont.)]			
Eastern (N.B., N.S., N.L.)	1.35	1.16 – 1.58	< 0.01*
Northern and Western (B.C., Alta., Sask., N.W.T.)	1.02	0.90 – 1.15	0.8
Charlson Comorbidity Index Value [Reference: 0]			
> 0	1.47	1.11 – 1.91	< 0.01*

Notes:

Excerpt of table. Only risk factors that were found to be statistically significant were included in this table.

* Statistically significant at $p < 0.05$.

‡ Other includes bearing surfaces that were reported as “other”, as well as smaller groups (for example, ceramic-on-metal and metal-on-ceramic).

§ Jurisdictions were combined into regional centres as follows: Western and Northern: British Columbia, Alberta, Saskatchewan and the Northwest Territories; Central: Manitoba and Ontario; Eastern: New Brunswick, Nova Scotia and Newfoundland and Labrador. This study does not include data from Quebec, Prince Edward Island, Yukon and Nunavut.

MoP: metal-on-polyethylene

MoM: metal-on-metal

Sources:

Canadian Joint Replacement Registry and Discharge Abstract Database, 2003–2004 to 2010–2011, CIHI

Conclusions

- Increased risk of revision for large-diameter metal-on-metal THAs and patients younger than 55.
 - Consistent with studies that have shown that the increased femoral head size in metal-on-metal implants may affect circulating metal ion levels.
 - May also suggest higher activity levels after primary procedure may contribute to early revision.
- Overall, findings from this study were consistent with those found in literature and by other international registries
- Largest study of its kind in Canada.
- Studies of implant attributes and other factors influencing revision rates important to reduce burden in health care system

Challenges Ahead

- Currently, lack of industry standardization in implant product labelling and barcoding
 - Difficult to accurately and systematically capture medical device product data
 - CJRR currently building product library and initiating work with stakeholders (manufacturers, GS1 Canada) to move towards more standardized labelling and product data capture.
- CJRR is the only national data source associating implant product with patient record.
 - Mandated in BC, MB, and ON (74% coverage nationally) but still voluntary in most regions
 - With increasing participation, more studies can be conducted to help improve the outcomes of hip and knee replacement patients



Thank you!



Canadian Institute
for Health Information
Institut canadien
d'information sur la santé